

Let's face it together!

The 5 most common mistakes made during wheelchair seating evaluations

Edward Milner



Background:

Occupational Therapist with a passion for Complex Rehab Technology (CRT) and neurological disability.

Experience:

Experienced in a wide range of clinical presentations from acute hospital in-patient management through to community rehabilitation services.

Proudly responsible for Education & Product Training around Australia for Medifab.

1. Skipping the MAT



© Confidential & Copyright Medifab

MAT Assessment



The Mechanical Assessment Tool (MAT) is a physical, musculoskeletal examination of the client and it includes:

1. Posture analysis on current seating system

- To help understanding the **problem**
- To determine current **loading areas** and **positioning strategies**


2. Supine assessment on the plinth

- To determine if deformities are **reducible** or not.
- To determine **angles** for sitting – hip flexion and knee extension are crucial!


3. Sitting assessment on the plinth – hand simulation

- To determine **optimal posture**
- To determine where **support surfaces** are needed

Translating information



Joint range of motion




Body segment angles



Seating angles

Assessment on the plinth

Client photo removed



Client photo removed

ROM of hips and knees

Pelvis-trunk- head
relationship

© Confidential & Copyright Medifab


Optimal (supported)
sitting posture

MAT Assessment


The MAT Assessment takes time and effort BUT it is critical:

- To determine available range for a seated posture;
- To determine seated angles;
- To determine loading areas;
- To determine whether postural asymmetries are reducible or non-reducible;
- To determine if neutral alignment is achieved, and what is the most optimal position for function
- To determine where support surfaces are required? Where? Of what kind?

Alternatives for community Ax.



© Confidential & Copyright Medifab



2. Jumping to product before defining parameters

Jumping to product before defining parameters

It can be tempting to start trying out equipment immediately, HOWEVER, you may end up:

- Addressing the problem but not the cause of it;
- Spending extra time as you didn't narrow down the options;
- Transferring the client multiple times to adjust the equipment or to try different solutions;
- Not being client/family centered.

Clinical reasoning worksheet

REFERRAL Client's current issues	OBJECTIVES Goals to achieve	CLIENT Problems and potentials		PRODUCT PARAMETERS Ways to achieve	POSSIBLE SOLUTIONS Products for trial
<p>1. Constant sliding which is impacting comfort, skin integrity and assistance care provision</p> <p>2. Reduced seating tolerance. E. can only sit for short periods of time.</p> <p>3. Bulky and heavy wheelchair which is currently impacting the way E. is transported on a vehicle</p>					

Clinical reasoning worksheet

REFERRAL Client's current issues	OBJECTIVES Goals to achieve	CLIENT Problems and potentials		PRODUCT PARAMETERS Ways to achieve	POSSIBLE SOLUTIONS Products for trial
<p>1. Constant sliding which is impacting comfort, skin integrity and assistance care provision</p> <p>2. Reduced seating tolerance. E. can only sit for short periods of time.</p> <p>3. Bulky and heavy wheelchair which is currently impacting the way E. is transported on a vehicle</p>	<p>✓ Sitting upright for long periods of time,</p> <p>✓ Sliding less on wheelchair,</p> <p>✓ Reduce manual handling,</p> <p>✓ Compact mobility base (for ease of transport and maneuverability)</p>				

Clinical reasoning worksheet

REFERRAL Client's current issues	OBJECTIVES Goals to achieve	CLIENT Problems and potentials	PRODUCT PARAMETERS Ways to achieve	POSSIBLE SOLUTIONS Products for trial
<p>1. Constant sliding which is impacting comfort, skin integrity and assistance care provision</p> <p>2. Reduced seating tolerance. E. can only sit for short periods of time.</p> <p>3. Bulky and heavy wheelchair which is currently impacting the way E. is transported on a vehicle</p>	<p>✓ Sitting upright for long periods of time,</p> <p>✓ Sliding less on wheelchair,</p> <p>✓ Reduce manual handling,</p> <p>✓ Compact mobility base (for ease of transport and maneuverability)</p>	<p>P E L V I S</p> <p>Non-reducible left pelvic obliquity.</p> <p>Slightly reducible right pelvic rotation.</p>		
		<p>H I P S</p> <p>90° hip flexion on the left; 80° hip flexion on the right.</p>		
		<p>L E , S</p> <p>Both knees tolerate 90° thigh-to-lower leg angle when available hip ROM is respected</p>		
		<p>T R U N K</p> <p>Hardly reducible scoliosis convex to left with pronounced rib hump.</p>		
		<p>H E A D</p> <p>Head control but strong tendency to rotate it towards the left.</p>		

Clinical reasoning worksheet

REFERRAL Client's current issues	OBJECTIVES Goals to achieve	CLIENT Problems and potentials	PRODUCT PARAMETERS Ways to achieve	POSSIBLE SOLUTIONS Products for trial
<p>1. Constant sliding which is impacting comfort, skin integrity and assistance care provision</p> <p>2. Reduced seating tolerance. E. can only sit for short periods of time.</p> <p>3. Bulky and heavy wheelchair which is currently impacting the way E. is transported on a vehicle</p>	<p>✓ Sitting upright for long periods of time,</p> <p>✓ Sliding less on wheelchair,</p> <p>✓ Reduce manual handling,</p> <p>✓ Compact mobility base (for ease of transport and maneuverability)</p>	<p>P E L V I S</p> <p>Non-reducible left pelvic obliquity.</p> <p>Slightly reducible right pelvic rotation.</p>	<p>Positioning cushion with the ability to accommodate left pelvic obliquity; also with the ability to accommodate 90° thigh-to-trunk angle on the left and 100° on the right, to accommodate for limited hip flexion.</p>	
		<p>H I P S</p> <p>90° hip flexion on the left; 80° hip flexion on the right.</p>	<p>4-point hip belt set up to accommodate hip flexion limitation (right) and pelvic rotation (left).</p>	
		<p>L E ' S</p> <p>Both knees tolerate 90° thigh-to-lower leg angle when available</p> <p>hip ROM is respected</p>	<p>Foot plates set up at 90°. Right one lower due to hip flexion limitation.</p>	
		<p>T R U N K</p> <p>Hardly reducible scoliosis convex to left with pronounced rib hump.</p>	<p>Deep contour back support that can accommodate scoliosis and rib hump.</p> <p>Lateral supports adjusted to an angle to support thoracic spine and promote optimal trunk alignment.</p>	
		<p>H E A D</p> <p>Head control but strong tendency to rotate it towards the left.</p>	<p>Posterior head support with lateral angle adjustment to correct head tendency.</p>	

Clinical reasoning worksheet

REFERRAL Client's current issues	OBJECTIVES Goals to achieve	CLIENT Problems and potentials	PRODUCT PARAMETERS Ways to achieve	POSSIBLE SOLUTIONS Products for trial
<p>1. Constant sliding which is impacting comfort, skin integrity and assistance care provision</p> <p>2. Reduced seating tolerance. E. can only sit for short periods of time.</p> <p>3. Bulky and heavy wheelchair which is currently impacting the way E. is transported on a vehicle</p>	<p>✓ Sitting upright for long periods of time,</p> <p>✓ Sliding less on wheelchair,</p> <p>✓ Reduce manual handling,</p> <p>✓ Compact mobility base (for ease of transport and maneuverability)</p>	<p>P E L V I S</p> <p>Non-reducible left pelvic obliquity.</p> <p>Slightly reducible right pelvic rotation.</p>	<p>Positioning cushion with the ability to accommodate left pelvic obliquity; also with the ability to accommodate 90° thigh-to-trunk angle on the left and 100° on the right, to accommodate for limited hip flexion.</p>	...
		<p>H I P S</p> <p>90° hip flexion on the left; 80° hip flexion on the right.</p>	<p>4-point hip belt set up to accommodate hip flexion limitation (right) and pelvic rotation (left).</p>	...
		<p>L E ' S</p> <p>Both knees tolerate 90° thigh-to-lower leg angle when available</p> <p>hip ROM is respected</p>	<p>Foot plates set up at 90°. Right one lower due to hip flexion limitation.</p>	...
		<p>T R U N K</p> <p>Hardly reducible scoliosis convex to left with pronounced rib hump.</p>	<p>Deep contour back support able to accommodate scoliosis and rib hump.</p> <p>Lateral supports adjusted in angle in order to support thoracic spine and promote optimal trunk alignment.</p>	...
		<p>H E A D</p> <p>Head control but strong tendency to rotate it towards the left.</p>	<p>Posterior head support with lateral angle adjustment to correct head tendency.</p>	...

3. Providing poor documentation



© Confidential & Copyright Medifab

Providing poor documentation

Not “painting a picture” can impact the success of your process.

Seating Assessment Form



spex

Client Data Assessment Form for Wheelchair Seating

Client Details:

Client Name: _____ Client Email Address: _____
 Date of Birth: _____ Care Name: _____
 Gender: _____ Care Contact Number: _____
 Client Contact Number: _____ Assessment Date: _____

Reason for Referral:

Medical History:

Diagnosis: _____
 Condition: ☐ Stable ☐ Deteriorating
 Cognition and perception: _____
 Medication: _____
 Hearing: ☐ Normal ☐ Impaired ☐ Deaf
 Vision: ☐ Normal ☐ Impaired ☐ Blind
 Respiration: ☐ Normal ☐ Ventilator dependent
 Sensation: ☐ Intact ☐ Impaired
 History of pressure injury (PI)? ☐ No ☐ Yes _____
 Risk of PI related to seating? ☐ No ☐ Yes _____
 Pressure relief: ☐ Independent ☐ Dependent ☐ Assisted
 Method: _____
 Pain history: _____
 Area of concern: _____
 Severity (1 = no pain, 5 = severe pain): _____

Postural Assessment:

Current wheelchair base: ☐ Stroller ☐ Manual ☐ Tilt-in-space
☐ Powerchair ☐ None
 Impact of current seating in ADLs: _____
 Postural needs and concerns: _____
 Transfers: ☐ Independent ☐ Needs assistance
☐ Dependent
 Client expectations from assessment: _____

Draw any asymmetries that need to be addressed in seating:

FRONT VIEW (FRONTAL PLANE) SIDE VIEW (SAGITTAL PLANE) TOP VIEW (TRANSVERSE PLANE)

Considerations for Outcome in Seating:

spex

Supine Assessment on the Plinth (MAT Evaluation)

Pelvis

Anterior Tilt: ☐ Reducible ☐ Non-Reducible
 Posterior Tilt: ☐ Reducible ☐ Non-Reducible
 Rotation: ☐ Left ☐ Right ☐ Reducible ☐ Non-Reducible
 Obliquity (Lateral Flexion): ☐ Left ☐ Right ☐ Reducible ☐ Non-Reducible

Hips

Dislocated/Subluxed: _____
 Flexion: _____
 Abduction: _____
 Adduction: _____
 Internal Rotation: _____
 External Rotation: _____

Knees

Flexion: _____
 Extension: _____

Feet

Dorsi-flexion: _____
 Plantar flexion: _____

Upper Limb

Shoulder Flexion: _____
 Elbow Flexion/extension: _____
 Wrist/hand: _____

Skin Inspection

Pelvis/buttocks: _____
 Trunk: _____
 Lower limbs: _____
 Upper limbs: _____

Muscle Tone

☐ Normal ☐ Increased ☐ Decreased ☐ Mixed
 Body segments affected: _____
 Describe: _____

Considerations for Outcome in Seating:

spex

Sitting Simulation on the Plinth

Sitting Balance

☐ Hands-free sitter ☐ Hands dependent sitter ☐ Proppped sitter

Pelvis

Neutral: ☐ Anterior pelvic tilt: ☐ Posterior pelvic tilt: ☐
 Rotation: ☐ Left ☐ Right
 Obliquity: ☐ Left ☐ Right

Trunk

Neutral: ☐ Scoliosis: ☐ Kyphosis: ☐ Lordosis
 Describe: _____

Head and Neck

Neutral: ☐ Extension: ☐ Flexion
 Describe: _____

Optimal position for sitting:

Seating Measurements:

LEFT		RIGHT	
A Lower Leg Length		F Chest Depth	
B Thigh Depth		G Hip Width	
C Ischial Well Length		H Chest Width	
D Shoulder Height		I Shoulder Width	
E Ankle Height		J Elbow Height	

Other Biomechanical Measurements:

Wheelchair Measurements:

LEFT		RIGHT	
B Frame Width		S Seat Depth	
T Back Support Cane Height		Footplate Hanger Angle	


Mobility Base Considered: ☐ Stroller ☐ Manual ☐ Tilt-in-space
☐ Powerchair ☐ None

Wheelchair Brand/Type: _____
 Wheelchair Cane Type: _____
 Wheelchair Notes: _____

Introducing the Shaperator!

A ground-breaking platform which assists in **prescribing**, **configuring**, **visualising**, and **ordering** Assistive Technology.

Instead of traditional ordering tools such as script forms and ordering guides, the Shaperator combines an easy-to-use ordering tool with 3D visualisation and intelligent product selection guides to take the pain out of the process and deliver the solution you're expecting.





<https://app.shaperator.com/login>

4. No follow-up

You did a great job so far! Don't stop there.

Follow up is essential to understand what is working versus what is NOT working!

- For client
- For support worker / caregiver

Work on strategies that work for everybody.

Use photos, sketches, videos, reminders, ...


5. Not considering Postural Care

We spend 1/3 of our lives in bed!

All the effort you have put into a new wheelchair process can unravel if the client is spending 8-10 hours a night in a destructive position.

Wheelchair seating must be part of an individualized 24-h Positioning Plan.

Posture Analysis Graph



Promote symmetry out of the wheelchair!

Well targeted interventions outside of the wheelchair can promote reduction of the postural tendencies seen in seating.

- Set up realistic goals. Can they be linked with seating?
- Promote **gentle** and **gradual** correction of destructive postural tendencies
- Aim for **supine lying** with pelvis leveled
- Support knees, trunk and feet accordingly.

Wheelchair Seating can be complex!

- Be **systematic**. Use forms, checklists and tools to support information collection and decision-making process.
- Be **pragmatic**. Discuss and set up realistic goals.
- Stay **informed**. Know what is available to support your clients better.
- Be **meticulous**. Provide in-depth and accurate information to other stakeholders.
- Promote **24-h Positioning**. Wheelchair seating is just part of it!
- Be patient and persevere!!

Thank you for attending!

Edward Milner, BSc OT

Clinical Educator | Medifab

ed.milner@medifab.com

The Medifab logo is located in the bottom right corner of the slide. It consists of the word "medifab" in a lowercase, sans-serif font, with a registered trademark symbol (®) to its upper right. The logo is white and is set against a dark blue circular background that has a white, torn-edge-like border on its left side.

medifab®